

## Development of a Standardised Analytical Method for the Detection and Quantitation of Melamine and Related Substances in Aqueous Media – Short Method Description –

Measurements are performed by high-performance liquid chromatography (HPLC) and MS/MS detection. The applied instrumentation and analytical conditions are listed in Table 3. Table 4 shows the resulting retention times and MS parameters for each target analyte and the available isotope-labelled internal standards.

Table 3: Summary of the analytical conditions.

HPLC system	HPLC Infinity 1290 (Agilent Technologies)
Injection volume	10 µL
Separation column	Luna HILIC column (150 x 2 mm, 3 µm, Phenomenex) with HILIC security guard cartridge (4 x 2 mm, Phenomenex)
Column temperature	30 °C
Eluents	A: acetonitrile B: ultrapure water + 10 mM ammonium acetate + 0.1% formic acid
Elution programme	Isocratic elution for 10 min: 90% A, 10% B
Flow rate	0.4 mL/min
MS/MS system	API5500 Triple Quad (Sciex)
Interface	Electrospray ionisation (ESI)
Ionisation mode	Positive (MEL, AMN, CYRO, HMMM) and negative (CYA, AMD)
Detection mode	Multiple reaction monitoring (MRM)
Mass transitions	See Table 4

Table 4: Retention times and MS parameters for detection of target analytes and internal standards.

Substance	Ionisation mode	Retention time (min)	Q1	Q3	DP (V)	EP (V)	CE (V)	CXP (V)
MEL	Positive	3.1	127	68* 85	106	10	37 25	10 12
AMN	Positive	3.5	128	86* 43.1	96	10	23 55	12 6
AMD	Negative	2.5	127	83.9* 42	-50	-10	-16 -28	-7 -5
CYA	Negative	1.5	127.9	42.1* 85	-45	-10	-26 -14	-5 -7
CYRO	Positive	2.0	167	68* 85	91	10	43 25	10 12
HMMM	Positive	1.2	391.1	177.1* 359.2	31	10	39 11	12 26
IS-MEL- <sup>13</sup> C3	Positive	3.1	130	87.1*	120	10	25	6
IS-AMN- <sup>13</sup> C3	Positive	3.5	131	88* 71	131	10	23 37	10 12
IS-AMD- <sup>13</sup> C3	Negative	2.4	130	86* 43	-55	-10	-16 -36	-5 -5
IS-CYA- <sup>13</sup> C3	Negative	1.5	130.9	43.1*	-45	-10	-38	-9
IS-CYRO-d4	Positive	2.0	171	86*	71	10	27	10

The method has been validated concerning recovery of evaporation, calibration linearity, matrix effects and reproducibility. Stability tests showed no analyte degradation in glass, PE, or PS containers over 28 days, regardless of the storage temperature. For drinking water, the method results in the limits of quantification (LOQ) given in Table 5. In other matrices, LOQ may vary depending on necessary sample dilution.

*Table 5:* LOQ values of target analytes in matrix drinking water.

<b>Substance</b>	<b>LOQ with direct injection</b>	<b>LOQ with evaporation</b>
MEL	0.10 µg/L	0.010 µg/L
AMN	0.20 µg/L	0.010 µg/L
AMD	0.10 µg/L	0.010 µg/L
CYA	1.0 µg/L	0.050 µg/L
CYRO	0.10 µg/L	0.010 µg/L
HMMM	0.10 µg/L	0.010 µg/L